

## How It Works



### THE SESSION

This simple, fast and effective treatment opens a new stay-young future. The treatment is unique, innovative, non-invasive, and let the pores absorb the nutrients and vitamins of our Coolifting serum better by the cold mist application

### APPLICATION ZONES

The combined cold CO2 gas flow and the proprietary hyaluronate facial solution with mild pressure can be applied around the eyes, forehead, cheeks, brow and zones around the lips.

It produces a more radiant facial appearance and an encourages skin hydration and attenuates fine wrinkle lines.



### NUMBER OF SESSIONS

The Coolifting protocol consists of:

- Single session: Immediate effect in a single session, a more radiant appearance and encouraged hydration of the skin.

It can be the beginning of a multiple session treatment with a more lasting effect.

- Multiple treatment protocol: 3-4 sessions (every 3 week) for a more lasting youthful look, increased radiant appearance, encouraged hydration of the skin.